




# LABORATORY LINES

LAMINAR FLOW CABINET



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



## KALSTEIN® Laminar Flow Cabinet

### Double person laminar flow cabinet

### Model Series YR05267



#### Product Description

Widely used in medical, college and other industries. provide air purification for the working area; is the essential equipment to carry out biotechnology research and experimental.

Mirror stainless steel, corrosion-resistant; low consumption, high efficiencies filter; microcomputer control, weight-balance front door.

Two people two side, two people single side; optional horizontal laminar flow and vertical laminar.

#### Standard Features

- ✓ Two models horizontal laminar flow and vertical laminar clean bench.
- ✓ According to using function are desktop and vertical clean bench, two people single side (two sides).
- ✓ The working surface is made of mirror stainless steel, easy to clean, corrosion-resistant.
- ✓ Equipped with many kinds filter to meet the customer's purification degree requirements.
- ✓ LED microcomputer controller, adjustable fan speed.
- ✓ Operation door can be located any position.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
 KALSTEIN FRANCE, S. A. S.



## KALSTEIN® Laminar Flow Cabinet

### Double person laminar flow cabinet

#### Model Series YR05267



#### Technical Parameters

Model	YR05267
Wind speed (m/s)	0.25 ~ 0.45
Illumination (LX) & Noise	Noise ≤62dB
Voltage (V)	220
Power (KW)	0.5
Inner size (mm) (H×W×D)	500×1250×430
Overoll size (mm) (H×W×D)	1610×1450×580
Packing size (mm) (H×W×D)	1200×1530×660
Net weight (KG)	120
Shipping weight (KG)	127



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed  
**The 3 articles of the week!**  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW  
**How does the climate affect your microscope?**  
microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through the use of lenses or a system of mirrors that collect and focus light from the object being observed.

**Kalstein Blog**  
**Article of the week:**  
How does the climate affect your microscope?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW  
**What is the importance of the laboratory refrigerator?**  
Refrigerator is a device that is one of the most common pieces of equipment that is used in laboratories. It is a controlled environment (refrigerated space), which allows you to store and maintain samples at low temperatures.

**Kalstein Blog**  
**Article of the week:**  
What is the importance of the laboratory refrigerator?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Allianc...**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

