




# LABORATORY SECTOR

REFRIGERATORS AND FREEZERS



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Refrigerators and Freezers**  
**316L 2-8°C Upright Medical Grade Refrigerator Pharmaceutical**  
**Model YR05288**



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



# KALSTEIN® Refrigerators and Freezers

## 316L 2-8°C Upright Medical Grade Refrigerator Pharmaceutical

### Model YR05288



#### Product Description:

YR05288 316L 2~8 C Upright Refrigerator, Pharmacy Refrigerator Medical Refrigerator for Plasma, Drugs, Vaccine etc.

#### Safety:

- ✓ Equipped with a complete temperature alarm system featuring audible buzzer and visible flashing light
- ✓ Capable of alerting failures due to high and low temperature sensor error, door ajar, , power failure and low battery

#### Humanized Design

- ✓ Optimized space utilization design with multi – level and adjustable shelf height to accommodate storage requirement of pharmaceutical products of various kinds. Safety lock to prevent unauthorized access
- ✓ Caster design and leveling feet
- ✓ LED Interior light

#### Reliability and Features

- ✓ High efficiency specialty compressor with known field reliability
- ✓ Permanently lubricated cooling fan for safety and longevity
- ✓ Forced- air cooling with optimized air distribution system designed
- ✓ To achieve maximum temperature uniformity and stability
- ✓ Optimized refrigeration system design for more effective cooling and
- ✓ Speedy recovery
- ✓ Temperature variation within  $\pm 3^{\circ}\text{C}$
- ✓ Adjustable temperature range 2-8°C
- ✓ Microprocessor control, digital display, and temperature adjustment with an increment of 0.1°C
- ✓ Large digital display screen



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



**KALSTEIN® Refrigerators and Freezers**  
**316L 2-8°C Upright Medical Grade Refrigerator Pharmaceutical**  
**Model YR05288**



**Technical Description:**

<b>Model</b>	YR05288
Performance	
<b>Temperature Range</b>	+2~+8°C
<b>Ambient Temperature</b>	10~32°C
<b>Controller</b>	Microprocessor
<b>Display</b>	Digital Display
Refrigeration System	
<b>Compressor</b>	1
<b>Refrigerant</b>	R290
<b>Cooling Type</b>	Forced-air cooling
<b>Defrost Type</b>	Frost free
Alarm	
<b>High/Low Temperature</b>	Y
<b>Power Failure</b>	Y
<b>Door ajar</b>	Y
<b>Thermostat failure</b>	Y
Accessories	



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**





# KALSTEIN® Refrigerators and Freezers

## 316L 2-8°C Upright Medical Grade Refrigerator Pharmaceutical

Model YR05288



<b>Caster</b>	4
<b>Lamp</b>	LED
<b>Door</b>	1/Glass door
<b>Shelf</b>	5
<b>Test hole</b>	1/25mm
<b>Temperature recorder</b>	Optional
Electrical Data	
<b>Power Supply(V/Hz)</b>	220/50-60; 115/60
<b>Rated Power(W)</b>	240W
<b>Electrical Current(A)</b>	1.2
<b>Power Consumption</b>	1.78kWh/24h
Dimensions	
Capacity	316L
Interior Size(W*D*H)	540*425*1380mm
Exterior Size(W*D*H)	660*555*1915mm
Package Size(W*D*H)	710*640*1970mm
N.t./G.t.(approx)	92.5/100kg
LoadingQty(20'/40'/40H)	27/54/54



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>




**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**1**  
Read this ARTICLE NOW



**How does the climate affect your microscope?**

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**2**  
Read this ARTICLE NOW



**What is the importance of the laboratory refrigerator?**

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Allianc...**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

