


Kalstein ® A different accompaniment, at your service

LABORATORY LINES ELECTROPHORESIS





All rights reserved ® KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,

• 34000 Montpellier, France.

Tiff +33 467158849 /+33 680760710/+33 663810

TIf: +33 467158849 / +33 680760710/ +33 66381002 https://kalstein.eu

KALSTEN FRANCE S. A. S.





OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.





#Letsgivemore 💗

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.





Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Basic Power Supply

Series YR03410



Basic Power Supply

YR series Basic power supply comprises 2 models, including: YR03410 and YR03411, which can meet the conventional horizontal electrophoresis, small and medium-sized vertical electrophoresis, membrane acetate electrophoresis experiment requirement, the lightweight and compact, simple and practical characteristics, make it popular in for basic teaching.

YR03410 Basic Power Supply



Type of Output: Constant voltage or constant current
Output range: JY300: 5~300V, 1~300mA, maximum 90W

Increment: 1V, 1mA

Timing Range: 1min~9hr, 59min

Display: LED

Output jacks: Two sets of output jacks
Dimension (L×W×H): 280×237×118 (mm)

Weight: 3.0 (kg)

- · Molding shell, touch keys, microprocessor intelligent control;
- · Automatic memory;
- · Standard timing operation;
- · Automatic detection of no-load, over-load.

YR03411 Basic Power Supplies



Type of Output: Constant voltage or constant current

Output range: JY600: 5~600V, 1~200mA, maximum 120W

Increment: 1V, 1mA

Timing Range: 1min~9hr, 59min

Display: LED

Output jacks: Two sets of output jacks
Dimension (LxWxH): 280x237x118 (mm)

Weight: 3.0 (kg)

- · Molding shell, touch keys, microprocessor intelligent control;
- · Automatic memory;
- Standard timing operation;
- Automatic detection of no-load, over-load.



#Letsgivemore 💗

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., Optimum Business Center 450 Rue Baden Powell,

34000 Montpellier, France.
 Tift: +33 467158849 / +33 680760710/ +33 663810023
 https://caldein.eu

KALSTEIN FRANCE, S.A. S





ANY DOUBT?Contact us!

Articles Kalstein

(Social Networks and Official Websites)







Article of the week: How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles





Article of the week: What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles

MONTPELLIER - FRANCE HEADQUARTER

Optimum Business Center, 450 Rue Baden Powell, 34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49 Tlf: +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/

Official Websites for Latin America.




Kalstein Training Videos

(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago