



# LABORATORY LINES

## MICROSCOPES



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.




## Binocular Infinity Plan Laboratory Biological Microscope Model YR05787



### Product Description

- ✓ Inclined 30°, Rotatable 360°, Interpupillary Distance 48-75mm, High Eyepoint Eyepiece,
- ✓ Which Is Convenient For People With Different Sitting Height; The Integrated Design Of The Eyepiece Lens Barrel,
- ✓ Avoids The Loss And Contamination Of The Eyepiece And Enhances The Integrated Performance Of The Microscope
- ✓ 4 Sets Objective Available: Achromatic, Plan, Infinity E-Plan, Infinity Plan.
- ✓ Mechanical Stage Size: 145 mm×140mm, Travel: 76mm×52mm, Scale: 0.1mm
- ✓ Abbe Condenser N.A.1.25, Iris Diaphragm, Rack & Pinion Up/Down
- ✓ 3W-LED Illumination Systems.



Abbe Condenser N.A.1.25, Iris Diaphragm, Rack & Pinion Up/Down



3W-LED Illumination Systems



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S.A. S



## Binocular Infinity Plan Laboratory Biological Microscope

### Model YR05787



Model	Cata. No.	YR05787	YR05788	
<b>Head</b>	Seidentopf Binocular Head Inclined 30°, Rotatable 360°, Interpupillary Distance 48-75mm	YR05788-op1	•	○
	Seidentopf Trinocular Head Inclined 30°, Rotatable 360°, Interpupillary Distance 48-75mm, Light Split Fixed P20:E80	YR05788-op2	○	•
<b>Nosepiece</b>	WF10X/20mm (High Eyepoint)	YR05788-op3	••	••
	Quadruple, Backward	YR05788-op4	•	•
	Quintuple, Backward	YR05788-op5	○	○
	PLAN 2.5x/0.07, W.D.=8.47mm	YR05788-op6	○	○
	PLAN 4X/0.10, W.D.=12.1mm	YR05788-op7	•	•
<b>Infinity Plan Objectives</b>	PLAN 10X/0.25, W.D.=4.64mm	YR05788-op8	•	•
	PLAN 20X/0.40(S), W.D.=2.41mm	YR05788-op9	○	○
	PLAN 40X/0.66(S), W.D.=0.65mm	YR05788-op10	•	•
	PLAN 60X/0.80(S), W.D.=0.33mm	YR05788-op11	○	○
	PLAN 100X/1.25(S, Oil), W.D.=0.12mm	YR05788-op12	•	•
	PLAN 100X/1.15(S, Water), W.D.=0.19mm	YR05788-op13	○	○
<b>Working Stage</b>	Mechanical Stage Size: 145 mm×140mm, Travel: 76mm×52mm, Scale: 0.1mm, Double Slide Holder	YR05788-op14	•	•
	Rackless (Integrated) Stage, Size 140×135mm, Travel Range 75×35mm, Double Slide Holder	YR05788-op15	○	○
<b>Focusing</b>	Coaxial Coarse & Fine Focusing, Coarse Focusing Travel Range: 26mm, Fine Focusing Travel Range: 2µm/Division	YR05788-op16	•	•
<b>Condenser</b>	Abbe Condenser N.A.1.25, Iris Diaphragm, Rack & Pinion Up/Down, With Filter Holder	YR05788-op17	•	•
	3W-LED Illumination Systems	YR05788-op18	•	•
<b>Illumination</b>	Powered By Safety Low Voltage Charger, Input Wide Voltage 100V-240V, Output 5V1A, Support Power Bank Supply For Outdoor Use	•	•	—
	12V/20W Halogen Lamp, Power Supply Wide Voltage 100V-240V	YR05788-op20	○	○
<b>Field Diaphragm</b>	6V/30W Halogen Lamp, Power Supply Wide Voltage 100V-240V	YR05788-op21	○	○
	Upgrade To Kohler Illumination	YR05788-op22	○	○
<b>Filter</b>	Dia.45mm, Blue	YR05788-op23	○	○
	Dia.45mm, Green	YR05788-op24	○	○
	Dia.45mm, Yellow	YR05788-op25	○	○
<b>Dark Field</b>	Dark Field Condenser, Dry, N.A.0.83-0.91, For 4x10x40x Objective	YR05788-op26	○	○
<b>Polarizing</b>	Analyzer + Polarizer Set	YR05788-op27	○	○
	Infinity Plan Phase Contrast Objective 10x + Phase Contrast Slide 10x + Phase Contrast Condenser + Telescope 11x + Green Filter	YR05788-op28	○	○
	Infinity Plan Phase Contrast Objective 20x + Phase Contrast Slide 20x + Phase Contrast Condenser + Telescope 11x + Green Filter	YR05788-op29	○	○
	Infinity Plan Phase Contrast Objective 40x + Phase Contrast Slide 40x + Phase Contrast Condenser + Telescope 11x + Green Filter	YR05788-op30	○	○
	Infinity Plan Phase Contrast Objective 100x + Phase Contrast Slide 100x + Phase Contrast Condenser + Telescope 11x + Green Filter	YR05788-op31	○	○
<b>Phase Contrast</b>	Turret Phase Contrast Attachment	YR05788-op32	○	○
	Infinity Plan Phase Contrast Objective 10x20x40x100x, Turret Phase Contrast Condenser 5-Holes, Telescope 11x, Green Filter	○	○	—
	C-Mount 0.5x, Focus Adjustable	YR05788-op33	○	○
	C-Mount 0.75x, Focus Adjustable	YR05788-op34	○	○
<b>Mirror Attachment</b>	Reflecting Mirror & Holder	YR05788-op35	○	○
<b>Package Weight</b>	Wire Winding Device	YR05788-op36	○	○
	Carton Size 390×280×450mm, 1 Unit/Ctn			
	Gross Weight: 7.5 kgs, Net Weight: 6.5 kgs			

Note: "•" In Table Is Standard outfits. "○" Is Optional Accessories



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación  
Jacinto Convit**  
<https://www.jacintoconvit.org/>


 **ONE  
TREE  
PLANTED**

<https://onetreepanted.org>

 **Humatem**  
Solidarité et coopération  
biomédicale  
<https://www.humatem.org/>

 **MANIAPURE  
FOUNDATION USA**  
<https://www.maniapurefoundation.org/>

**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

refrigerator is a device that is one of the most common pieces of laboratory equipment. Its purpose is to maintain, in a controlled environment (refrigerated space), various fluids and substances

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?


Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago

