

LABORATORY SECTOR

REFRIGERATORS AND FREEZERS







All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S.A.S





OUR SERVICES Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings In any part of the world, receive your

induction or training from our specialized team of engineers



Quick Response Our work team is always available to response all your consults or questions, in order to support you in any situation.





#Letsgivemore 💙

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





Kalstein HARD'S TRANSFER Ξ 三三三日 ΞΞΞ EEEE = = EEEEE Ξ = = BBBBB $\equiv \equiv \equiv$ Ξ $\Xi \equiv \Xi$ Ξ EHH 1 11 11 11 a 3 6) 61



#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S CE

All rights reserved © KALSTEIN France S. A. S.,



Product Description:

656L Double Door Pharmacy and Lab Refrigerator with Glass Door and LED Interior Light

Ratability and Key Features:

- ✓ High efficiency specialty compressor with known field reliability
- ✓ Permanently lubricated cooling fan for safety and longevity
- ✓ Forced- air cooling with optimized air distribution system designed
- ✓ To achieve maximum temperature uniformity and stability
- ✓ Optimized refrigeration system design for more effective cooling and

Safety:

- ✓ Speedy recovery
- ✓ Temperature variation within ±3°C
- ✓ Adjustable temperature range 2-8°C
- ✓ Microprocessor control, digital display, and temperature adjustment with an increment of 0.1°C
- ✓ Large digital display screen

Ergonomic Design:

- ✓ Optimized space utilization design with multi level and adjustable shelf height to accommodate storage requirement of pharmaceutical products of various kinds.
- ✓ Safety lock to prevent unauthorized access
- ✓ Caster design and leveling legs
- ✓ LED Interior light



#Letsgivemore Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.

All rights reserved © KALSTEIN France S. A. S.,



Equipped with a complete temperature alarm system featurin

✓ Capable of alerting failures due to high and low temperature

sensor error, door ajar, , power failure and low battery

g audible buzzer and visible flashing light



All rights reserved ® KALSTEIN France S. A. S. 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



www.kalstein.eu



Technical Description:

Performance		
Temperature Range	2~8°C	
Ambient Temperature	+10~+32°C	
Controller	Microprocessor	
Display	Digital Display	
Refrigeration System		
Compressor	1	
Refrigerant	R600a	
Cooling Type	Forced air cooling	
Defrost Type	Automatic	
Alarm		
High/Low Temperature	Y	
Sensor Error	Y	
Power Failure	Y	
Door ajar	Y	



#Letsgivemore Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



www.kalstein.eu



Accessories		
Caster	4	
Light	LED	
Shelf	10	
Temp. Recorder	Qu'in a l	
Chart recorder	Optional	
Test hole	1/25mm	
Electrical Data		
Power Supply(V/Hz)	220/50~60; 115/60	
Rated Power(W)	385	
Electrical Current(A)	2.3	
Power Consumption	5.75 kWh/24h	
	Dimensions	
Capacity	656L	
Interior Size(W*	D*H) 1100*454*1325mm	
Exterior Size(W	• D • H) 1240*630*1895mm	
Package Size(W	*D*H) 1335*710*2000mm	
N.t./G.t.(appro	>>>) 166/184kg	
Loading Qty(20'/40	5'/40'H) 13/27/27	



#Letsgivemore Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.

All rights reserved © KALSTEIN France S. A. S.,



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



www.kalstein.eu





#LETSGIVEMORE

WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org





ANY DOUBT? Contact us!



Articles Kalstein (Social Networks and Official Websites)





PARIS - FRANCE **HEADQUARTER** 2 Rue Jean Lantier 75001 Paris - France

Tif: +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/

Kalstein Training Videos (Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes 445 views • 1 week ago

