




MEDICAL LINES

ELECTROCARDIOGRAPH



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics


We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





Electrocardiograph

Model YR7031

- ✓ ECG wave form 30000 patient save/replay
- ✓ 7 inch color & touch screen
- ✓ Lead:standard 12 leads
- ✓ 3 channel ECG machine with interpretation
- ✓ Manual/Auto/Arrhythmia Analysis Modes
- ✓ 80mm 3-ch recording with rhythm lead
- ✓ Built in rechargeable Li-ion battery
- ✓ Adaptor 110-230V, 50/60Hz AC power supply
- ✓ Built in USB/RS232 interface(ECGNET software Optional)
- ✓ YR7031 (3 channel)/ YR7062 (6 channel)

Power Source	Electric
Warranty	Lifetime
After-sale Service	Return and Replacement
Material	metal, plastic
Shelf Life	1years
Quality Certification	ce
Instrument classification	Class II
Safety standard	None
Type	Cheap Price Animal ECG
Display	7 color LCD Touch Screen Super-slim design
Sale points	Support SD card, USB disk
Analysis Modes	Manual/Auto/Arrhythmia
Lead	12 leads
ECGNET software	Optional
Anti-baseline Drift	Automatic
Interface	USB/RS232
Storage	30000 patient save/replay
Recording	112mm *20m

Patient cable	1set
Limb electrode	4sets
Chest electrodes	6sets
Grounding cable	1set
Paper axis	1set
Thermal paper	1set
Fuse	2sets
3-cord power cord	1set
AC Adaptor	1set



#Let'sgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.




All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Titl: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Electrocardiograph

Model YR7031

Features	7 "display, touch screen Support Ecg analysis capabilities
Lead	standard 12 leads
lead acquisition	synchronously 12 leads
Input circuit	Floating; Protection circuit against Defibrillator effect
Input Impedance	$\geq 50\text{M}\Omega$
Input circuit current	$\leq 0.05\mu\text{A}$
Record mode	Automatic: 3CH \times 4+1R, 3CH \times 4, 3CH \times 2+2CH \times 3,3CH \times 2+2CH \times 3+1R,6CH \times 2; Manual: 3CH, 2CH, 3CH+1R, 2CH+1R; Rhythm: Any lead selectable.
Filter	EMG Filter: 25 Hz / 30 Hz / 40Hz/75 Hz / 100 Hz / 150Hz DFT Filter: 0.05 Hz/ 0.15 Hz AC Filter: 50 Hz / 60Hz
CMRR	>100dB
Patient current leakage	<10 μA
Input Circuit Current	<0.05 μA
Frequency Response	0.05Hz~150Hz (-3dB)
Sensitivity	2.5mm/mV, 5 mm/mV, 10 mm/mV, 20 mm/mV
Anti-baseline Drift	Automatic
Time constant	$\geq 3.2\text{s}$
Noise level	<15 $\mu\text{Vp-p}$
Paper speed	12.5 mm/s, 25 mm/s, 50 mm/s
Recording mode	Thermal printing system



#Letsgivemore


Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form an image of the object being observed.
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator in a laboratory is one of the most important pieces of equipment. It is used to maintain, in a controlled environment (refrigerated space), various fluids and substances.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



MONTEPELLIER - FRANCE
SEDE

Optimum Business Center,
450 Rue Baden Powell,
34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49

Tel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

Official Websites for Latin America.

