




# LABORATORY LINES

## ANALYZER



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Analyzer**  
**Hematology Analyzer For Veterinary**  
**Model YR05123-V**



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



# KALSTEIN® Analyzer

## Hematology Analyzer For Veterinary

### Model YR05123-V



#### Features:

- ✓ Animal Type: cat, dog, horse, mouse, rat, rabbit, pig, cow, monkey, sheep, user-defined settings
- ✓ 3-part differentiation of WBC, 24 parameters, double chamber with 60 sample tests per hour
- ✓ Advanced valve technology, long life
- ✓ Electrical resistance for counting and SFT method of Hemoglobin
- ✓ 8.4 color touch screen, Windows interface all testing parameter displayed simultaneously
- ✓ Automatic diluting, Mixing, rinsing and clog clearing
- ✓ Automatically sample probe cleaning (inside and outside)
- ✓ Large storage capacity: more than 100,000 samples +3 histograms
- ✓ Internal! thermal-sensitive printer or external! printer.

#### Parameters

- ✓ 3 histograms: WBC,RBC and PLT
- ✓ RDW-SD, PLT, MPV,PDW,PCT,L-PCR,P-LCC
- ✓ WBC, Lymph#, Mid#, Gran#, Lymph%, Mid%, Gran%, RBC, HGB,HCT, MCV, MCH, MCHC, RDW-CV

#### Specifications:

- ✓ Methodology: Electrical resistance for counting, hemoglobin cyanide method and SFT method for hemoglobin
- ✓ Parameter: 3-part differentiation of WBC; 21 parameters and 3 color histograms (WBC, RBC, PLT) Work mode: Double Channel + unique Hemoglobin test system
- ✓ Sample volume: 9.6µL for Venous and capillary mode, 20µL for Prediluted mode Throughput: More than 60 samples per hour, auto sleeping and waking-up functions
- ✓ Storage: More than 100,000 sample results including histograms can be stored, convenient for inquiry and management of history data
- ✓ Operation System: Windows operation system graphical buttons mouse and keyboard operation is available Operation language: English, Spanish, French, Portuguese and Indonesian
- ✓ QC control: X-B, L-J, X, SO, CV %
- ✓ Alarm: 25 test of alert messages
- ✓ Input/ output : RS232, parallel printer and keyboard
- ✓ Print: Graphic thermal printer with various printing format, optional external! printer Temperature: 18°( - 30°( ,wet :5 10-90%
- ✓ Power supply: 220 V ±22 VAC, 50±1Hz Dimension: 33 CM ( L ) \* 38 CM ( W ) \* 43 CM ( H ) Weight: 20 KG

#### Precision:

Parameters	linear range	CV%
WBC( $10^9 / L$ )	0.0 - 99.9	:52%
RBC ( $10^{12} / L$ )	0.0 - 9.99	:51.5%
MCV ( fl )		:50.5%
PLT ( $10^9 / L$ )	40-150	:54.0%
HGB( g/L )	0-999	:51.5%
	0.0-300.0	



#### #Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación  
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>




**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**1**  
Read this ARTICLE NOW



**How does the climate affect your microscope?**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**2**  
Read this ARTICLE NOW



**What is the importance of the laboratory refrigerator?**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Allianc...**

23 views • 3 years ago



**Kalstein's Microscopes**

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

