


# Kalstein A different accompaniment, at your service

## **LABORATORY SECTOR**

**PIPETTES** 





All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S





**OUR SERVICES** 

### **Benefits and Support**

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



## Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



#### **Quick Response**

Our work team is always available to response all your consults or questions, in order to support you in any situation.





#### #Letsgivemore \*\*

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



#### **Technical Support**

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.





#### **Delivery Logistics**

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



#### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



## **KALSTEIN® Pipettes**

### **Liquid Handling Pipettes**

Model YR03732 // YR03732-8







#### #Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S



#### **KALSTEIN®** Pipettes

#### **Liquid Handling Pipettes**

Model YR03732 // YR03732-8






#### **Product Description:**

YR Pipette Series are designed and manufactured to provide you with a range of robust, accurate and precise pipettes. These models are a fixed volume, air-displacement pipette with the selected volume shown on a digital indicator (volume display)

#### Features:

- ✓ 3 digits
- ✓ Magnifying shape
- ✓ Chemical resistance
- ✓ 121°C Sterilization
- ✓ Very low operating force
- ✓ Color indicates pipette volume
- ✓ Positioned for perfect ergonomics
- Only a few turns to reach your desired volume



Ejector

Very low operating force
 Positioned for perfect ergonomics

Model	Volume Range	Accuracy	Repeatability
YR03732	0.1 – 2µl	±12.0 – 1.5 %	≤6.00 – 0.70%
YR03732-1	0.5 – 1µl	±2.5 – 1.0 %	≤1.50 − 0.40%
YR03732-2	2 – 20µl	±2.5 – 1.0 %	≤1.50 – 0.30%
YR03732-3	5 – 50μΙ	±2.0 – 1.0 %	≤1.50 – 0.25%
YR03732-4	10 – 100μl	±1.8 – 0.8%	≤0.50 – 0.15%
YR03732-5	20 – 200µl	±1.8 – 0.6%	≤0.50 – 0.15%
YR03732-6	100 – 1 000µl	±1.5 – 0.6%	≤0.30 – 0.15%
YR03732-7	500 – 5 000μl	±1.2 – 0.6 %	≤0.30 – 0.15%
YR03732-8	1 – 10ml	±2.5 – 0.6 %	≤0.30 – 0.15%



#### #Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S. 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S







#LETSGIVEMORE >

## WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org



Solidarité et coopération biomédicale

https://www.humatem.org/



https://www.maniapurefoundation.org/



## **ANY DOUBT?**Contact us!

#### **Articles Kalstein**

(Social Networks and Official Websites)







Article of the week: How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles





Article of the week: What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles



PARIS - FRANCE HEADQUARTER

2 Rue Jean Lantier 75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/

#### **Kalstein Training Videos**

(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago

