




LABORATORY LINES

MICROPLATE READER
AND WASHER



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Microplate Washer

Model YR05124 & YR05125



Features:

- ✓ 8-well and 12-well wash head
- ✓ One or more strips can be chosen by keyboard directly
- ✓ Compatible with Flat, U-or V-bottom plates
- ✓ Shaking or soaking, time can be adjustable
- ✓ Waste and lotion and distilled water liquid warning automatically.

Specifications:

- ✓ Cv <2%
- ✓ Washing head:
- ✓ 8- and 12-well strips, Interchain
- ✓ Residual Volume:
- ✓ Residual Volume: 1µL/well
- ✓ Certificates: CE ISO9001 ISO13485
- ✓ Washing programs:
- ✓ 100 programs can be stored
- ✓ Washing Cycles:
- ✓ 1-99, optional
- ✓ Soaking Time:
- ✓ 0-999s, adjustable
- ✓ Dispense volume:
- ✓ 50µl-950µl, adjustable



Remark:

The model YR05125 has a built-in internal-temperature cyclic ELISA Plate pallet incubator (patent). Other's specifications are same as the model YR05124.



Dispense volume	0-950uL, adjustable
Washing Cycles	1-99 optional
Washing Mode	1-12, row and rank, selectable
Soaking time	0-999s adjustable
Vibration time	0-999s adjustable
Display	Large LCD
Channel	2 channels(including one channel for water)
Memory Interface	UP to 100 programs can be stored RS232
Power supply	220V±10% 50/60Hz AC110V±10% 50/60Hz 100W
Dimension	440mm(L)×305mm(W)×154mm(H)
Net Weight	13KG
Micro-plate Types	48/96-well U-/V-/Flat-bottom plates or strips



#Letsgivingmore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





Microplate Washer Model YR05126

Features:

- ✓ Centrifugal dripping, no clapping.
- ✓ No clotting, without crossing contamination
- ✓ Two plate washing at the same time.

Specifications:

- Washing head: 96 needles
- Washing wells: 24, 48, 72, 96 wells optional, can continuous arrange the wells, no interval.
- Micro-plate type: Flat, Round, U, V
- Washing mode Single plate (A or B) double plate (A-B)
- Lotion: A, B
- Washing time : 1-999 times adjustable, step 1s
- Flushing time : 0-999 seconds adjustable, Step s1
- Pressure setting : 0-9.9 adjustable
- Centrifugal time : 0-999 seconds adjustable, Step 1s
- Centrifugal rate : 1200-2500rpm
- Residual liquor: 0.1 uL per well
- Washing head position: 0- 10.00 mm continuously adjustable, completely washing
- Liquid level induction : insufficient washing and distilled Water alarm
- Capacity bottle : 2pcs 10000mL bottles for lotion, 1 pc for distilled water, another for waste
- 7 inch color touch screen operation, convenience and efficient with low failure rate.
- Power Supply: SMPS: 110/220V ± 10%, 50/60Hz ±1HZ
- Working Condition: Temp. 5 C° -40 C°; Relative humidity ≤ 80 %
- Dimensions : 51cm (L)*51cm(W)*43cm
- Weight : 30KG



- ✓ Single needle recoil, waste out of the plate automatically, no blocking. Clean Plate thoroughly, not waste plate.



- ✓ Centrifugal dripping, no clapping. Transparent cover, the cover safety inspection.



- ✓ Washing two plates for one time. Washing one plate cost less than 25 seconds.



- ✓ Tinting the base, automatic liquid discharge, no waste in the bucket



- ✓ Vertical plate washing method, put into the plate upward with 60°, to avoid the outflow



- ✓ 7 inch touch screen, digital adjustable washing time and pressure.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tit: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Microplate Reader

Model YR05127

Specifications

- Microplate types: 96-well plate(U,V, or flat-bottom)
- Wavelength: 405nm,450nm,492nm,630nm,4 free positions for additional filters
- wavelength accuracy < 2nm
- Reading range: 0.000–4.000 Abs
- Measurement range: 0.000–3.000Abs
- Light source: LED, >100,000 hours
- Measuring system: 8-channel optical system
- Test mode: End-point, kinetic, Multi-wave analysis
- Resolution <1% or ± 0.001 Abs
- Stability: $< \pm 0.005$ Abs/10min
- Well Coincidence: ± 0.015 A
- Linearity: $r > 0.999$
- Repeatability: < 0.3%
- Reading speed: 5s single wavelength , 10s double wavelength
- QC: West guard Multi-rule, automated QC alert.
- Self-Checking function: Self-checking and calibration function.
- CPU: High-speed CPU
- Printer: External jet printer, multi-form patient report.
- Interface: RS232 interface, USB.
- Power supply: 220VAC 15% or 110VAC 15%50–60Hz 100W
- Work environment: Temp 5 \pm 40 \pm ; Humi \pm 85%;
- Machine Dimension: 44CM*37CM*24CM
- Weight: 13 KG



Features:

High performance Portable Elisa Microplate Reader 96 well plate.

- ✓ Windows operation interface, figure button, mouse operation, touch screen.
- ✓ 96-well plate, multiple tests in the single plate.
- ✓ Arrangement of Test item in micro-plate can be store and re-use.
- ✓ Parameter can be modified and stored.
- ✓ Plate shaking, speed and time adjustable.
- ✓ Open system, test item and parameters can be added or modified.
- ✓ Large storage for 1000 test item and 1,000,000 results.
- ✓ 10.4 inch colour touch LCD.
- ✓ Printing comprehensive report, include hospital and patient information and test results.
- ✓ Built-in portable computer, easy to use.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Microplate Reader

Model YR05128



Features:

Low cost but high performance multi functional Automated Elisa Plate Reader.

- ✓ Windows operation interface, figure button, mouse operation.
- ✓ 96-well plate, multiple tests in the single plate.
- ✓ Arrangement of test item in micro-plate can be store and re-use.
- ✓ Plate shaking, speed and time adjustable.
- ✓ Open system, test item and parameters can be added or modified.
- ✓ Large storage for 1000 test item and 1,000,000 results.
- ✓ External large LCD colour screen.
- ✓ Printing comprehensive report, include hospital and patient information and test results.
- ✓ External computer, easy to use.



Specifications:

Microplate types: 96-well plate(U, V, or flat-bottom)
Test mode : End-point, Kinetic, Multi-wave analysis
Measuring system:8-channel optical system
Light source:Tungsten Halogen, 12V/20W,>2000 hours
Wavelength: 405、450、492、630nm, 4 more filters optional.
wavelength accuracy : < 2nm
Reading range:0.000-4.000 Abs
Measurement range: 0-3.00Abs
Resolution : 0.001Abs (display) , 0.0001Abs(Inner counting)
Stability: <±0.005Abs/10min
Sensitivity: ≥ 0.01L /mg
Linearity: r > 0.999
Repeatability: < 0.3%
Reading speed : 5s single wavelength , 10s double wavelength
QC:Westguard Multi-rule,automated QC alert.
Self-Checking function : Self-checking and calibration function.
Printer: External jet printer, multi-form patient report.
CPU: High-speed CPU
Interface :RS232 interface,USB, PSTN and PC connection.
Power supply: 220VAC±15% or 110VAC±15% 50-60Hz 100w
Work environment: Temp 5℃—40℃ : Humi ≤85%
Machine dimension : 43CM*34CM*16CM
Weight : 7 KG



#Let'sgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tel: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form...

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator is a device that is one of the most common pieces of equipment used in laboratories to maintain, in a controlled environment (refrigerated space), various fluids and substances...

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



MONTPELLIER - FRANCE
HEADQUARTER
Optimum Business Center,
450 Rue Baden Powell,
34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49
Tif: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

Official Websites for Latin America.

