




LABORATORY LINES

BALANCES



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Balances Multifunction Weighing Scale



Function:

- ✓ Appearance design and utility model access to patents
- ✓ 7 inch touch screen
- ✓ Large range and high precision
- ✓ Large size stainless steel pan
- ✓ AC and DC power supply
- ✓ Date/Time set
- ✓ Tare function/Counting
- ✓ Check weigher/Percentage/Density
- ✓ Net and Gross/Sum
- ✓ Unit conversion(g/ct/oz)...
- ✓ Overload alarm
- ✓ Connect PC and printer
- ✓ Optional label printer
- ✓ Function Expansion:Under weighing/
- ✓ Remote display/Internal calibration (YR05615 & YR05617 option)



connect to computer/printer/USB



Percentage

Convenient for customers to read the formula ratio

Density

Can test solid/liquid
Need to be equipped with density frame



USB flash disk programming function

Thermal printer form

Check weighing upper and lower limit setting

It is suitable for production enterprises to check and weigh finished products to improve work efficiency



#Let'sgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Balances

Multifunction Weighing Scale



Specifications					
Model	YR05615	YR05616	YR05617	YR05618	YR05619
Capacity	10kg	30kg	10kg	20kg	30kg
Readability	1g		0.1g		
Operating Temp	5°C-35°C				
Repeatability	±3d				
Linearity	±3d				
Stable Time	≤2S				
Sensitivity Drift	±3d				
Corner	±3d				
Weighing Units	g、oz、ct...				
Pan Size	305x230mm				
Key Name	Touch Screen				
Appearance Size	395x118x320mm				
Packing Size	460x365x210mm				
Gross Weight	7.5kg				
Display	7Inch Touch Screen				
Interface	RS232/RJ45/USB...				
Internal Calibration	option	/	option	/	/
Weighing Mode	Counting/Check Weighing/Percentage/Density				
Function Expansion	Under Weighing/Remote Display				
Power	Input:220V 50HZ/110V 60HZ Output:9V				
Battery	6V/4AH Rechargeable Battery				



#Let'sgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 


**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>


 **ONE
TREE
PLANTED**

<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>

 **MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects at a very small scale. This is achieved through a complex system of lenses and mirrors, which form...

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator is a device that is one of the most common pieces of equipment used in laboratories. It is a controlled environment (refrigerated space), various fluids and substances...

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

