




MEDICAL SECTOR

OPERATION TABLES



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Operation Tables
Operation Tables
Model YR02153



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S.



KALSTEIN® Operation Lamps**Operation Tables****Model YR02153****Product Description:**

Impressive improvements on functional performance, user interaction and operational safety make YR02153 a truly universal operating table, which offers effective solutions for the widest variety of requirements in the OR. It is an ideal operating table that can be used in general surgery, thoracic surgery, obstetrics and gynecology, orthopedics, urology, brain surgery, and ENT surgery.

Features:

- ✓ 350mm Longitudinal Shift: Adding more space for C-arm in X-ray imaging surgeries, such as cardio surgery, orthopedic, ERCP.
- ✓ Double Back Plates: Kidney elevator structure can be reached by double back plates, with 0-105mm adjustable height.
- ✓ Stable Valley-shaped Base: Solid and stable valley-shaped base is suitable for low-position MIS surgeries, like plastic, neurosurgery and ophthalmic surgeries.
- ✓ Alternative Pedal Control: Providing more convenience and efficiency especially in laparoscopy and endoscopic surgery.
- ✓ Override Control: Backup control panel on the column offers more safety in power supply interruption.
- ✓ 500mm Tabletop Height Adjustment Range: this product is universally applicable, and especially suits for Neurosurgery, ENT in the 600mm low-height-position.
- ✓ Intuitive Handle Control: Adjusting height and tabletop positions, including the back plate using the backlit keypad.

**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Operation Lamps**Operation Tables****Model YR02153****Technical Description:**

Model	YR02153
Table length	2140mm
Table width	520mm
Table height	1100mm / 600mm
Trendelenburg / Reverse	30 ° / 30 °
Lateral tilt (right / left)	25/25 °
Head plate (up / down)	45 ° / 90 °
Back plate (up / down)	70 ° / 40 °
Leg plate (up / down / split)	10 ° / 90 ° / 90 ° x 2
Flexion / Reflection	110 ° / 1220 °
Operation load	200KG
Input power	AC 110 / 220V, 60 / 50Hz

**#Let'sgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

