




MEDICAL LINES

OPERATING TABLES



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





YR02158

Mechanical-hydraulic Operating Table

Kalstein's YR02158 model is a versatile operating table, adopting safe & stable hydraulic drive and pneumatic-assisted system, realizing all kinds of surgical positioning required clinically, with simple manual control.



Kidney & Thoracic Surgery

MIS (Minimal Invasive Surgery)

Urologic & Gynecological Surgery



Neurosurgery



General Surgery



Anorectal Surgery



Four-stand base of great stability. High mobility with front-ended release and brake control.



Integrated kidney elevator 0-120mm height adjustment.



High-quality hydraulic system is adopted, which can realize required positioning, like tilting and lifting.

YR02158 Specification

Table Length	2140mm
Table Width	500mm(not include side rails)
Table Height	1040mm/690mm
Trendelenburg/Reverse	20°/15°
Lateral Tilt (left/right)	20°/20°
Head Plate (up/down)	30°/90°
Back Plate (up/down)	75°/20°
Leg Plate (down/split)	90°/90° x 2
Flex/Reflex	192°/105°
Operation Load	200KG



#Let's Give More

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S., Optimum Business Center 450 Rue Baden Powell, • 34000 Montpellier, France.
 Tit: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Select Guide



Functionality	Purpose & Application	Model
Longitudinal Shift	Adding more space for in C-arm X-ray imaging surgeries, such as cardiosurgery, orthopedic, ERCP, lithotomy.	YR02153, YR02154
Interchangeable Head and Leg Section	Interchangeable leg and head section give bigger irradiated area for C-arm use.	YR02153, YR02155, YR02156, YR02157
Low Height Position	Wide tabletop height adjustment range with low height position (550~600mm) is suitable in neurosurgery, ENT and plastic surgeries.	YR02153, YR02157
Double Back Sections	Kidney elevator structure can be reached by double back plates, used in kidney surgery.	YR02153
Kidney Elevator	Build-in or external kidney elevator make it easier for positioning in kidney surgery.	Build-in: YR02153, YR02155/6/7, YR02158 External: YR02154
Pedal Control	Providing more convenience and efficiency especially in laparoscopy and endoscopic surgery.	YR02153, YR02155
Override Control	Backup control panel on the column offers more safety in emergency.	YR02153, YR02154
Backlit Remote Handle	Making it easy to control the OT in dark OR ambient environment, such as during the endoscopic and laparoscopic surgeries.	YR02153
Motorized Back Plates	The back plates are electro-driven, suitable for lithotomy, neurosurgery, and other sitting position surgeries.	YR02153, YR02155
Valley-shaped Base	New design table base provides excellent stability and safety for surgeries, like neurosurgery, plastic surgery.	YR02153



#Letsgivemore


Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tit: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form...

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator is a device that is one of the most common pieces of equipment used in laboratories to maintain, in a controlled environment (refrigerated space), various fluids and substances...

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



MONTPELLIER - FRANCE
HEADQUARTER
Optimum Business Center,
450 Rue Baden Powell,
34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49
Tif: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

Official Websites for Latin America.

