



LABORATORY LINES

CENTRIFUGES



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



PRP Centrifuge Model YR469



Features;

- ✓ Microcomputer control, DC brushless motor, stable operation, low noise and high-speed accuracy.
- ✓ LED or LCD display. user-friendly interface, simple and convenient operations.
- ✓ Special program is made for PRP kit.
- ✓ LED or LCD display. user-friendly interface, simple and convenient operations.
- ✓ Various of PRP kit is available.
- ✓ Special program is made for PRP kit.
- ✓ The real-time conversation Of Speed and RCF is convenient for operation.

Technical performance

- ✓ Microcomputer control, DC brushless motor, stable operation, low noise and high-speed accuracy.
- ✓ It has self-locking device, lid safety device, over-speed safety device, and automatic alarm device, etc.
- ✓ Silastic sealed ring comfort to GMP.US FDA certified.

The PRP operation Demo



Blood collection
Transfer blood to
PRP Kit



Put PRP kit into the rotor



From the bottom to the top:
Red blood cell, PRP and plasma



PRP extraction



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tit: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



PRP Centrifuge Model YR469

CE EWC
ISO9001
ISO13485



FDA



Specification	YR469	YR469-1
Max Speed	5000r/min	4000r/min
Max RCF	4100xg	2220xg
Max Capacity	50mlx4	50mlx4
Speed Accuracy	± 30r/min	±30r/min
Time Setting	1min~99min	1min~99min
Noise	<65dB(A)	<65dB(A)
Power Supply	AC220V 50Hz 10A	AC220V 50Hz 5A
Dimension(LxWxH)	410mmx490mmx310mm	330mmx420mmx280mm
Packing Size(LxWxH)	500mmx580mmx400mm	430mmx520mmx390mm
Net Weight	22kg	14kg

Rotor



#Letsgivemore


Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form...

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator in a laboratory is one of the most important pieces of equipment. It is used to maintain, in a controlled environment (refrigerated space), various fluids and substances...

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago




MONTPELLIER - FRANCE
HEADQUARTER
Optimum Business Center,
450 Rue Baden Powell,
34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49
Tif: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

Official Websites for Latin America.





LABORATORY LINES

CENTRIFUGES



Kalstein®

A different accompaniment, at your service

