



# LABORATORY SECTOR

PH METERS



## Kalstein®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® PH Meters**  
**Pocket pH Tester**  
**Model YR01804**



**#Letsgivingore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



## KALSTEIN® PH Meter Pocket pH Tester Model YR01804



### Product Description:

Digital pocket pH meter, Manual Temperature Compensation. The tester comes with a plastic body pH electrode for outdoor applications, accuracy: 0.01pH.

### Features:

- ✓ BNC connector is capable of connecting the different types of pH electrodes.
- ✓ 1 to 3 points calibration with auto-buffer recognition.
- ✓ Manual temperature compensation allows adjusting the temperature value.
- ✓ Auto-Read function senses and locks the measurement endpoint.
- ✓ Auto-Power Off effectively conserves battery life.
- ✓ Setup menu allows setting the pH buffer group, number of calibration points, temperature
- ✓ unit, etc.
- ✓ Reset function automatically resumes all settings back to the factory defaults.

### Applications:

- ✓ Water and wastewater treatment
- ✓ Aquaculture and hydroponics
- ✓ Environmental monitoring and other industrial applications



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



**KALSTEIN® PH Meter**  
**Pocket pH Tester**  
**Model YR01804**



**Technical Description:**

<b>Model</b>	YR01804
<b>Range</b>	-1.00 to 15.00pH
<b>Resolution</b>	0.01pH
<b>Accuracy</b>	±0.01pH
<b>Calibration Points</b>	1 to 3 points
<b>pH Buffer Options</b>	USA (pH4.01/7.00/10.01), NIST (pH4.01/6.86/9.18)
<b>Temperature Compensation</b>	0 to 100°C/32 to 212°F, manual
<b>Hold Function</b>	Manual or auto-endpoint
<b>Auto-Power Off</b>	8 minutes after last key pressed
<b>Connector</b>	BNC
<b>Display</b>	Dual-line LCD (21×21mm)
<b>Power</b>	2×1.5V AAA batteries
<b>Battery Life</b>	Approximately 200 hours of continuous use
<b>Dimensions</b>	175(L)×40(Dia.)mm
<b>Weight</b>	100g



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed  
**The 3 articles of the week!**  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW  
**How does the climate affect your microscope?**  
Microscopes are instruments used for observing microscopic objects. This is a sensitive device that requires a controlled environment to maintain its performance.  
**Kalstein Blog**  
**Article of the week:**  
How does the climate affect your microscope?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW  
**What is the importance of the laboratory refrigerator?**  
Refrigeration is a key element of the most advanced equipment that is used in laboratories. It is a controlled environment (refrigerated space), without fluctuations in temperature.  
**Kalstein Blog**  
**Article of the week:**  
What is the importance of the laboratory refrigerator?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Alliances**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

