



# LABORATORY SECTOR

MICROSCOPE



## Kalstein®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



## KALSTEIN® Microscope Polarizing Microscope Model YR0261



### #Letsgivemore


Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



## KALSTEIN® Microscope Polarizing Microscope Model YR0261



### Product Description:

- ✓ Quintuple Nosepiece Center Adjustable
- ✓ Strain-free Infinity Plan Objectives  
4x10x40x100x
- ✓ Round Stage Dia.170mm Center Adjustable
- ✓ Analyzer Unit With Bertrand Lens,  $\lambda$  Slip,  $1/4\lambda$  Slip, Quartz Wedge
- ✓ Transmit Halogen 24V100W Brightness Adjustable



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



## KALSTEIN® Microscope Polarizing Microscope Model YR0261



### Technical Specifications:

<b>Model</b>	YR0261	
<b>Objective</b>	Strainfree Infinity Plan Achromatic Objective 20x/0.4/∞/- W.D. 5.1mm	YR0261-op1
<b>Stage</b>	Attachable Stage, Moving Range 30*40mm	YR0261-op1
<b>Filter</b>	Netural Filter ND25	YR0261-op1
	Netural Filter ND26	YR0261-op1
<b>Adapter</b>	Photo Adapter	YR0261-op1
	C-Mount 1.0x	YR0261-op1
<b>Items</b>	YR0261	
<b>Optical System</b>	Infinity Optical System	
<b>Head</b>	Seidentopf Trinocular Head, 30° Inclined, Interpupillary 48-75mm, Diopter Adjustable	
<b>Eyepiece</b>	Plan Eyepiece EW10x/22mm, Eyepiece Tube Dia.30mm	
	Plan Eyepiece With Scale of Crosshair	
	Plan Eyepiece With Cross	
	Plan Eyepiece With Gridding	
<b>Objective</b>	Strainfree Infinity Plan Achromatic Objective	
	4x/0.10/∞/- W.D.17.3mm	
	10x/0.25/∞/0.17 W.D.10mm	
	40x/0.65/∞/0.17 W.D.0.54mm	
<b>Nosepiece</b>	100x/1.25/∞/0.17 W.D.0.13mm	
	Backward Quintuple Nosepiece, Center Adjustable	
	360° Dial Rotation Analyzer, Minimum Scale Reading 0.1° (Vernial Scale)	
	Switch Between Orthoscopic And Conoscopic Observation	
<b>Analyzer</b>	Bertrand Lens Postion Adjustable	
	Optical Compensator: λ Slip (First Class Red), 1/4λ Slip, Quartz Wedge	
	Polarizing Rotatable Stage, Dia.170mm, Center Adjustable	
	360° Scale, Minimum Division 1°	
<b>Working Stage</b>	Minimum Reading 6' By Means of Vernial Scale	
	45° Click Stop Knob	
	Strainfree Swing Condenser N.A.0.95/0.25	
<b>Condenser</b>		
<b>Focusing</b>	Coaxial Coarse & Fine Adjustment, Focusing Range 32mm, Fine Division 0.001mm, Coarse Stroke 37.7mm Per Rotation, Fine Stroke 0.1mm Per Rotation	
<b>Polarizer</b>	360° Rotatable, 0° Position Adjustable	
<b>Filter</b>	Blue Filter, Dia.45mm	
<b>Illumination</b>	Transmit Illumination Halogen Lamp 24V100W, Brightness Adjustable	



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale


<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Allianc...**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

