

LABORATORY SECTOR

SPECTROPHOTOMETER







All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S.A.S





OUR SERVICES Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings In any part of the world, receive your

induction or training from our specialized team of engineers



Quick Response Our work team is always available to response all your consults or questions, in order to support you in any situation.





#Letsgivemore 💙

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





KALSTEIN® Spectrophotometer UV-VIS Spectrophotometer With Ultraviolet Visible Scanning Software Model YR01857





#Letsgivemore 💚

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved (® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



All rights reserved © KALSTEIN France S. A. S.,

www.kalstein.eu

KALSTEIN® Spectrophotometer UV-VIS Spectrophotometer With Ultraviolet Visible Scanning Software Model YR01857



Product Description:

YR01857 spectrophotometer equipped with 6 inches LCD display, is an ideal and advanced analytical instrument for laboratory to realize wavelength scanning, Kinetics test, multi wavelength functions. All functions can be operated on spectrophotometer and can be read directly on the display.

6 Inches Display

YR01857 series has a 6 inches LCD display to show results and curves directly on the screen.

8mm Thick Optical Base:

YR01857 uses a rigid 8mm die-case aluminum base as its optical mount to ensure instrument stability and reliability.

Powerful Functions:

Multi functions like spectrum scanning, standard curve, kinetics, multi wavelength, DNA/RNA/Protein testing can be operated directly on the spectrophotometer and all corresponding curves and data can be displayed directly.

Data Output

YR01857 is equipped with USB port to connect with a PC to display spectrum scanning, kinetics and multi wavelength testing results on the screen. The software is optional.

Model	YR01857		
Wavelength Range	190-1100nm		
Bandwidth	1.8nm		
Wavelength Accuracy	±0.5nm		
Wavelength Repeatability	≤0.2nm		
Photometric Accuracy	±0.3%T		
Photometric Repeatability	≤0.15%T		
Photometric Range	-0.3-3A, 0-200%T, 0-9999C		
Stability	± 0.002A/h @ 500nm		
Baseline Flatness	± 0.002A/h		
Noise	± 0.0005A		
Stray Light	≤0.05%T @ 220nm, 360nm		
Data Output Port	USB		
Printer Port	Parallel Port		
Display	320*240 Dots LCD		
Lamps	Tungsten Lamp & Deuterium Lamp		
Detector	Silicon Photodiode		
Power Requirements	AC 220V/50Hz or AC 110V/60Hz		
Dimension	460*380*180mm		
Weight	20kg		



#Letsgivemore 💗

All rights reserved © KALSTEIN France S. A. S.,

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S. 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



KALSTEIN® Spectrophotometer UV-VIS Spectrophotometer With Ultraviolet Visible Scanning Software Model YR01857



ltem	Description	Quantity	Unit
1	Spectrophotometer	1	set
2	1cm Glass cuvette	4	pcs
3	Power cord	1	pcs
4	User's Manual	1	pcs
5	1 cm quartz cuvette	2	pcs
6	Dust Cover	1	pcs



#Letsgivemore Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.

All rights reserved © KALSTEIN France S. A. S.,



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



www.kalstein.eu





#LETSGIVEMORE

WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org





ANY DOUBT? Contact us!



Articles Kalstein (Social Networks and Official Websites)





PARIS - FRANCE **HEADQUARTER** 2 Rue Jean Lantier 75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/

Kalstein Training Videos (Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes 445 views • 1 week ago

