

## LABORATORY LINES

**REFRIGERATORS AND** FREEZERS







All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S.A.S





# our services Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



#### Online Inductions and Trainings In any part of the world, receive your

induction or training from our specialized team of engineers



#### **Quick Response** Our work team is always available to response all your consults or questions, in order to support you in any situation.

n



### #Letsgivemore 💚

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



## Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



## **Delivery Logistics**

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.







Series YR05297 (Cascade System)









#Letsgivemore Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities. This product is specially designed and manufactured for long term storage of various biological product, including viruses, germs, erythrocytes, leucocytes and cutis. Applications can be found in blood banks, hospitals, epidemic prevention services, and research institutes, laboratories in electronic and chemical plants. Biological engineering institutes and marine fishery companies

#### Safety

- ✓ Malfunction alarms including high and low temperature, power failure, sensor error, clean-filter, and extremely high ambient, abnormal voltage, thermostat failure, low battery, condenser clean, door ajar.
- Capable of producing two types of alarm outputs, audible buzzer and visible flashing light.
- ✓ Door open feature standard and USB port for temperature date downloading standard on upright models.
- ✓ Remote alarm contacts.

#### **Installation and Application**

- ✓ Wide rage operating voltage system from 187V to 242V designed to allow units installed in areas with poor voltage condition.
- ✓ Suitable for 10°C to 32°Cambient temperature.

#### **Filed Proven Reliability**

- ✓ Unique insulated inner door design for four separate storage compartments to minimize frost buildup inside the chamber.
- ✓ Specialized control system design for a well-balanced operation of cascade refrigeration system.
- ✓ Positive filed proven reliability record.

#### **Key Design Features**

- ✓ Microprocessor-controlled system designed for controlled rage of -40°C to -86°C for cabinet space with 1°C increment
- ✓ Settable high temperature and low temperature alarms.
- ✓ Automatic clean-filter alarm and sensor error alert.
- ✓ Adjustable storage shelf height
- ✓ Optional temperature recorder, storage racks and storage box

#### **Energy saving**:

- ✓ Unique door seal design for the minimum loss of cold temperature during a door opening.
- ✓ High performance VIP insulation panel to minimize cabinet heat gain and to improve temperature stability.
- ✓ HC Refrigerant with better efficiency and less energy cost, good for the environment.



All rights reserved ® KALSTEIN France S. A. S. 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



www.kalstein.eu

LET'S

GIVE

MORES

## -40 °C To -86° Ultra Low Freezer 588L Series YR05297 (Cascade System)



CE



2 individual inner doors can be opened independently to minimize frost buildup inside the chamber Unique door seal design for the minimum loss of cold temperature during a door opening Stainless Steel handle to ensure the door open conveniently even in the case of frost



Excellent Doors Seals 5 gasket seal, 4 seals for out door

7 Inches touch Screen All data display on the screen, easy to observe

Write Board for short records which would pass to next user

Green Handle, easy to use with 2 locks

Special V.I.P(Vacuum Insulation Panel) design with extra insulation fatorcuts down heat loss to cabinet by 25%

-1) =

Big condenser for better heat dissipation

Changeable sensor box or you could put another sensor from here to monitor or record temperature





#Letsgivemore V Thanks to your purchase, a donation

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S. 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



All rights reserved © KALSTEIN France S. A. S.,

www.kalstein.eu

Series YR05297 (Cascade System)

#### 7 Inches Touch Screen & Functions






Energy-saving and enviromental friendly

CE









Alarm

High Temperature

Low Temperature

Power Failure

Probe Failure

Low Battery

Hot Condenser

Other alarms

High Ambient Temperature

Door Ajar

#### #Letsgivemore Thanks to your purchase, a donation is made

to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S

6



www.kalstein.eu


All rights reserved © KALSTEIN France S. A. S.,

Series YR05297 (Cascade System)





#### **Dual System Circulation System**



#Letsgivemore 💗 Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S.A.S



All rights reserved © KALSTEIN France S. A. S.,

www.kalstein.eu



Series YR05297 (Cascade System)



	Model	YR05297	YR05298	YR05299	YR05300	YR05301
Specificailons	Shelf/drawer QTY	Stainless steel shell/3	Stainless steel shel 1/3	Stainless steel shell/3	Stainless steel shell/3	Stainless steel shelf/3
	Cooling	Direct cooling	Direct cooling	Direct cooling	Direct cooling	Direct cooling
	Defrost	Manual	Manual	Manual	Manual	Manual
	Refrigerant	HC	HC	HC	HC	HC
	Power Consumption(kWh/24h)	12.8	11.6	12	15.4	15
	Noise (db)	50	50	53	53	53
	Ambient temperature (*C)	10~32°C	10~32°C	10~32°C	10~32°C	10~32°C
	Temperature range (*C)	-40~-86°C	-40~-86°C	-40~-86°C	-40~-86°C	-40~-86°C
Cooling system	Compressor brand/QTY	SECOP/2	SECOP/2	SEC OP/2	SEC OP/2	SECOP/2
	Sensor	PTC	PTC	PTC	PTC	PTC
	Temperature controller	Microprocessor	Microprocessor	Microprocessor	Microprocessor	Microprocessor
	Display	LCD display	LCD display	LCD display	LCD display	LCD display
Power	Voltage/frequency (V/Hz)	220/50Hz	220/50Hz	220/50Hz	220/50Hz	220/50Hz
	Power (W)	870W	845W	815W	935W	890W
	Current (A)	6.75A	7A	7A	7A	7A
Material	Inside material	Stainless steel	Stainless steel	Stainless steel	Stainless steel	Stainless steel
	Outside material	Color sprayed steel	Colorsprayed steel	Color sprayed steel	Color sprayed steel	Color sprayed steel
	Insulation	PURF	PURF	PURF	PURF	PURF
	Capacity (L/cu.ft)	338/11.93	408/14.40	588/20.76	728/25.70	838/29.59
Size	NT/GT. (kg)	250/305	292/317	307/368	360/401.5	397/430
	Exterior size(W*D*H)(mm)	880*980*1940	900*940*1940	880*1080*2000	1260*890*2010	1190*1045*2010
	Interior size(W*D*H)(mm)	488*607*1140	606*575*1180	606*738*1310	963*558*1350	894×718×1310
				975*1140*2200	1320*970*2190	1270*1105*2220
	Package size(W*D*H)(mm)	940*1070*2130	955*960*2100			
	CBM	2.14	1.93	2.45	2.80	3.12
	20GP/40GP/40HQ	12/24/24	12/24/24	12/24/24	8/18/18	8/18/18
Frozen Racks	Racks	12(4X4)	16(5X4,4X4)	16(5X5)	18(5x4)+6(6x4)	24(5X5)
	Boxes (2 Inches)	192	288	400	504	600
	Samples(2ml)	19200	28800	40000	50400	60000
Alarm	Hgh/low temperature	Y	Y	Y	Y	Y
	Powerfailure	Y	Y	Y	Y	Y
	Abnormal voltage	Y	Y	Y	Y	Y
	Sensor error	Y	Y	Y	Y	Y
	Filter screen check	Y	Y	Y	Y	Y
	Thermostat failure	Y	Y	Y	Y	Y
	Low battery	Y	Y	Y	Y	Y
	Condenser clean	Y	Y	Y	Y	Y
	High ambient temperature	Y	Y	Y	Y	Y
	Door ajar	Y	Y	Y	Y	Y
	Power failure backup		72h	72h	72h	72h
Accessories	system (a arm)	Y	Y	Y	Y	Y
	Inside door	2	2	2	2	2
	Test hole QTY/diameter	2/25MM	1/25MM	1/25MM	1/25MM	1/25MM
	USBport	Y	Y	Y	Y	Y
	Chart recorder	Optional	Optional	Optional	Optional	Optional
	Data logger	Optional	Optional	Optional	Optional	Optional
	Remote alarm port	Y	Y	Y	Y	Y
Optional	CO2 back up system	Optional	Optional	Optional	Optional	Optional
	LN2 back up system	Optional	Optional	Optional	Optional	Optional
Other	Certification	ISO.CE	ISO.CE	ISO.CE	ISO,CE	ISO,CE



#### #Letsgivemore 💗

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S

All rights reserved © KALSTEIN France S. A. S.,





#LETSGIVEMORE

## WITH THE **ADCQUISITION OF A KALSTEIN EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org

Humatem Solidarité et coopération biomédicale https://www.humatem.org/



https://www.maniapurefoundation.org/



## ANY DOUBT? Contact us!



(Social Networks and Official Websites)





PARIS - FRANCE **H E A D Q U A R T E R** 2 Rue Jean Lantier 75001 Paris - France

**Tif:** +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu


https://kalstein.eu/ https://www.kalstein.fr/

#### Kalstein Training Videos (Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes 445 views • 1 week ago

