


Kalstein A different accompaniment, at your service

LABORATORY LINES VORTEX MIXER







Tlf: +33 467158849 / +33 680760710/ +33 66381002 https://kalstein.eu

KALSTEN FRANCE S.A.S.





OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.





#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.





Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





Vortex Mixer

YR05761

Description

is best solution for shaking 96 well micro plate and other small size testing tube, centrifuge tube and pcr strip, etc. Applications in cell and molecular biology, microbiology and immunology, etc.

Feature

- 1. Intelligent controller for speed and timer;
- 2. Adopt permanent magnet rotor, with over-temperature protection, ensure long time working;
- 3. Brushless motor, noiseless and maintenance-free;
- 4.Two operation mode "touching" and "continuous", easy to switch;
- 5. Can automaticly identify various clamps quickly and ensure safe using;



Input power	AC100-250V/50/60Hz
Output power	18W
Drive mode	DC24V brushless motor
Control mode	PWM steplessly variable frequency
Operation mode	"jog", "continuous"
Power fuse	PPTC Resettable Fuses /no replacement
Maximum load	1.5Kg
Speed Range	0/200-3000r/min.
Oscillation Mode	Circumference
Circle diameter	4.5mm
Fixture Identification function	Automatic identification/ Safe and efficient
Speed function /showing	Full digital / bright LED display
Timing function /showing	Speed timing function /showing Digital timing 1S-999min
	countdown "jog" timing / bright LED display
Dimension(LxWxH)	198×138×67mm
Weight	4.1kg
The noise for Maximum speed	≤43dB(A)
Maximum acceleration	≤3s
Maximum deceleration	≤3s
Ambient temperature / relative humidity	5-40℃/80%
Protection grade	IP21
Testing standard	EN 61010-1,EN 61010-2-020,EN 61326-1,EN 61010-3-2/A2



#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved (t) KALSTEIN France S. A. S., Optimum Business Center 450 Rue Baden Powell,

34000 Montpellier, France.

Tif: +33 467158849 / +33 680760710/ +33 663810023

https://kaistein.eu Kalstein France, S.A. S



CE

Vortex Mixer YR05761



Optional working head



M8-1 standard

D<50mm tubes

Max speed: 3000rpm



M8-3 standard

ELISA plate

Max speed: 1500rpm



M8-2 for florence flask

D<85mm flask

Max speed: 3000rpm



M8-4 for testing tube

30x10mm

Max speed: 1500rpm



for testing tube

5x17mm

Max speed: 1500rpm



M8-6 for testing tube

16x12-13mm

Max speed: 1500rpm



M8-7 for testing tube

12x15mm

Max speed: 1500rpm



ਅਲ-ਲ for testing tube

0.2x8ml pcr strip or 96well microplate

Max speed: 1500rpm



M8-9 for testing tube

12x17mm

Max speed: 1500rpm



M8-10 for testing tube

2x50ml

Max speed: 1500rpm



M8-11 for testing tube

4x25mm penicillin bottle

Max speed: 1500rpm



#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved @ KALSTEIN France S. A. S., Optimum Business Center 450 Rue Baden Powell,

34000 Montpellier, France.

Tif: +33 467158849 / +33 680760710/ +33 663810023

https://kaistein.eu Kalstein France, S.A. S







#LETSGIVEMORE >

WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org



Solidarité et coopération biomédicale

https://www.humatem.org/



https://www.maniapurefoundation.org/




ANY DOUBT? Contact us!

Articles Kalstein

(Social Networks and Official Websites)







Article of the week: How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites





Article of the week: What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles

Kalstein Training Videos

(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



MONTPELLIER - FRANCE HEADQUARTER

Optimum Business Center, 450 Rue Baden Powell, 34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49 **Tlf:** +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/

Official Websites for Latin America.

