




MEDICAL LINES

WHEELCHAIRS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Electric Climb Stairs Wheelchair

Model YR05446



Characteristics:

- ✓ Long-range.
- ✓ Lithium wheelchair for stair-climbing wheelchair overall.
- ✓ It comes with a big lithium battery and is three times faster than other power wheelchairs.
- ✓ It also works well on all kinds of stairs regardless of the texture.

Description

Kalstein Electric climb stairs wheelchair YR05455, turns dreams into reality; with this unique device, you can explore the world like never before and bring a new taste to your everyday routines! The YR05455 wheelchair can help you climb stairs, overcome obstacles, and travel in any direction by using the one joystick on the side of the wheelchair. They are extremely popular because they're easy to use and spare the users from having to use their arms to propel the chair. It also makes them less dependent on another person. These chairs are meant for everyone and can be especially useful for people with more serious disabilities who can't propel a manual chair.

The newest designs come with some amazing features such as a speech device, various controls, options to adapt the body position, and more.



- **Electric Power move on the ground**
- **Electric power climb stair.**



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
 KALSTEIN FRANCE, S. A. S



Electric Climb Stairs Wheelchair

Model YR05446



Technical Description

Model	YR05446
Frame	Aluminum alloy
Battery	18Ah lithium battery
Wheel size	8 & 12 inches
Net weight	38kg
Max loading	150kg
Climb/Descent	three-speed
Climbing ability	500 stair steps
Stair width required	> 700mm
Stair slope required	<45 °
Stair depth required	1100mm
Product size	98*65*125cm
Folding size	110*65*35cm
Seat width	43cm
Seat depth	45cm
Cushion	Breathable sponge
Color	Black/Yellow/Red



#Let'sgive**more**

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



Manual Wheelchair

Model YR05455



Characteristics:

- ✓ PU integral wheel
- ✓ Comfortable honeycomb, breathable cushion
- ✓ Portable handbag, foldable & light for easy travel
- ✓ Four brake designs are smoothly advanced
- ✓ Carbon steel raw materials reinforcement design

Product Description

Kalstein medical manual Wheelchair is a series of cost-effective products. According to the needs of different people, there are many functions, such as turning over, lifting the back, lifting the legs, and releasing the intestines. At the same time, to meet the needs of hospital beds, we also launched a series of nursing beds for families and individuals, providing customers with more choices.



- **Backside brake system**
- **Removable footrest**
- **Manual break system**
- **Lifting armrest**



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



Manual Wheelchair

Model YR05455



Technical Description

Model	YR05455
Seat Width	42 cm
Seat Depth	43 cm
Seat Height	45 cm
Standing Height	140 cm
Backrest Height	40 cm
Overall Length	106 cm
Overall Width	67 cm
Overall Height	85 cm
Rear Wheel Size	24 inches (62 cm)
Front Wheel Size	5 inches (12 cm)
Net Weight	25.8 Kg
Gross Weight	33.6 Kg
Frame Material	Steel
Rear Wheel Material	Aluminium Alloy
Load	100 Kg
Package Size	115 x 65 x 64 cm



#Let'sgive**more**

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



Manual Wheelchair

Model YR05454



 **Kalstein**



Characteristics:

- ✓ Chromed frame
- ✓ Fixed foot rest
- ✓ TPR front wheel
- ✓ Pneumatic rear wheel
- ✓ Handle with brake

Product Description

Kalstein medical wheelchair is our most advantageous product. We can provide wheelchairs with different functions for people with different needs, including electric wheelchairs, manual wheelchairs, wheelchairs for children with cerebral palsy, standing wheelchairs, stair climbing wheelchairs and other wheelchairs with different functions.

- **Removable footrest**
- **Manual break system**
- **Steel Frame**

Technical Description

Model	YR05454
Seat Width	42 cm
Seat Depth	43 cm
Seat Height	45 cm
Standing Hight	140 cm
Backrest Height	40 cm
Overall Length	106 cm
Overall Width	67 cm
Overall Height	85 cm
Rear Wheel Size	24 inches (62 cm)
Front Wheel Size	5 inches (12 cm)
Net Weight	25.8 Kg
Gross Weight	33. 6 Kg
Frame Material	Steel
Real Wheel Material	Aluminium Alloy
Load	100 Kg
Package Size	115 x 65 x 64 cm



#Let'sgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



Manual Wheelchair

Model YR05453



Characteristics:

- ✓ Chromed frame
- ✓ Fixed armrest
- ✓ Fixed foot rest
- ✓ TPR front wheel
- ✓ Pneumatic rear wheel
- ✓ Drop back handle
- ✓ Portable handbag
- ✓ Handle with brake

Product Description

Kalstein medical manual Wheelchair is a series of cost-effective products. According to the needs of different people, there are many functions, such as turning over, lifting the back, lifting the legs, and releasing the intestines. At the same time, to meet the needs of hospital beds, we also launched a series of nursing beds for families and individuals, providing customers with more choices.

- **Backside brake system**
- **Come with commode**
- **Manual break system**
- **Footrest can be rubber or steel**



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
 KALSTEIN FRANCE, S. A. S



Manual Wheelchair

Model YR05453



Technical Description

Modelo	YR05453
Cushion	Heavy Duty Oxford Woven Cushion
Rear tire	Solid tire 24 inch alloy hub
Front tire	7 inch non-pneumatic
Width (open)	680mm
Width (fold)	260mm
Seat width	46mm
Length from floor to seat	500mm
Height for chair	900mm
Length for chair	1 000mm
Not tire	410mm
Size for package	94x26x96
Cushion	Heavy Duty Oxford Woven Cushion
rear tire	Solid tire 24 inch alloy hub
Front tire	7 inch non-pneumatic
Width (open)	680mm
Width (fold)	260mm



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
 KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

